

# BRO as a passenger seat.

This way you can use BRO as a passenger seat in the car.<sup>1</sup>

BRO has a special «ride-on mode»<sup>2</sup> that turns BRO into a passenger seat in the car. This guide only gives you an overview of how to use the passenger-mode in a wheelchair accessible vehicle. For more detailed instructions and other car models, please be sure to read the user manual. You can download the latest version from our website.

Still not sure how to use BRO as a passenger seat? We're happy to help you. You can find the contact details of us and our partners on our website.

1

Drive your BRO facing forwards in front of the loading ramp in order for you to enter through the boot towards the designated seating location within your vehicle and switch to parking mode. Now fasten the hooks of the safety belts<sup>3</sup> into the two front fastening points of your BRO or ask an assisting person to do this.



2

Now select the car hook icon on the touchpad to activate the passenger-mode.



3

Now switch to passenger-mode by selecting the track mode icon and drive forwards over the loading ramp into the boot.<sup>4</sup>



<sup>1</sup> Only use BRO in vehicles approved for transporting a wheelchair.

<sup>2</sup> Only wheelchairs from Series One (V1.1) are equipped with the attachment points.


<sup>3</sup> For your security, use belts that are approved for a weight of at least 160 kg (according to ISO 10542-1). We recommend securing it with 6 safety straps (2 front and 4 rear).

<sup>4</sup> On very steep ramps, ask an assisting person to make sure that BRO cannot tip over.

4

When you have reached the designated position in the car, select the parking mode icon to switch to the passenger position. If you move the joystick, your BRO will now automatically move to the recommended passenger position (ring lights up green).

If it is not possible to take the recommended passenger position due to the conditions of the vehicle, you can end the transition early. To do this, select the car hook icon on the touchpad. You can now adjust your seat position as usual (ring lights up yellow). We recommend that you stay in the pre-programmed position whenever possible.



5

Fasten yourself with a three-point harness<sup>1</sup> and ask an assisting person to secure your BRO at the two rear attachment points.



6

Once you have reached your destination, release the rear seat belts or ask an assisting person to do so and switch back to passenger-mode by selecting the track mode icon to drive out.


Be sure to leave the front safety belts fastened while driving out and ask an assisting person to make sure that BRO cannot tip over!



7

As soon as you have completely left the loading ramp, select the car hook icon on the touchpad to end the passenger-mode. Now remove the front safety belts.

Now you can switch to driving mode via the parking mode and continue your adventure.



<sup>1</sup>Never use the security belt of BRO here, but always the integrated belts of the car. Often an additional certified security belt is required as an extension.