



«Thanks to BRO, I can discover the world with my children independently and spontaneously.»

José



«I have a disability, but when I sit in BRO, it doesn't show.»

Anne





«With BRO, I can go about my daily work independently again and even long distances no longer matter.»

Jeanette


Book a non-binding test drive now and experience freedom and independence.

www.scewo.com/en/testdrive



 SWISS DESIGN. SWISS ENGINEERING. SWISS MADE.

Scewo AG | Technoparkstrasse 7 | 8406 Winterthur | Switzerland
info@scewo.ch | +41 44 500 86 03 | www.scewo.com | Q2 2023



scewo BRO

Scewo BRO App
Intuitive operation
Remote control mode
Additional settings



Requires **no upper-body stability**



35 km
range

Thresholds

Effortlessly overcoming thresholds with BRO.

Service Hotline



Anti-tip system for more safety when driving.

Up to 10 km/h
*depending on regulations



With two large wheels you can navigate **safely and nimbly** on various surfaces.

DRIVING MODE



1050 MM

688 MM

Joystick adapter for various special attachments.



72%
Track mode

Stairs up to 72 % are no problem for BRO.

10%
Driving mode

BRO effortlessly handles **gradients** up to 10 %.


162 KG

BRO incl. battery

TRANSPORT MODE




Partner network for competent advice and service near you.



Suitable and certified for transport in cars and aeroplanes.



Adapt BRO to your needs. Adjust the backrest and footrest at the touch of a button or make various mechanical adjustments to the headrest, armrest, footrest and more.



Stand-up aid & buttocks relief



Interchangeable seat and back cushions (also from third-party suppliers)



BRO is recognised by insurance companies.



Dynamic seat tilt




Thanks to **over-the-air updates** BRO is always up to date.



Add-ons

Thigh supports
Swivelling interface
Convenient car ramp

Multiple award-winning design



All technical data for BRO can be found on our website.



Sensors

detect the end of the stairs for automated stair climbing.



A rear view camera for greater safety when navigating.



Sahara Gold
Urban Grey
Polar White

